


MICRO SPRING NEEDLE HOLDER WITHOUT LOCK

17.5cm Straight | SKU 1379



PRODUCT OVERVIEW

Summary	Straight. Without lock. Overall length...
Category	Ophthalmic Instruments, Needle Hol...
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Castroviejo Length: 17.5 cm Instrument Type: Needle Holder Blade Type: Straight Tip Type: Blunt Handle: Thumb (double spring) Material: Stainless Steel Alias Names: Castroviejo microsurgical needle holder Surgical Instrument: Hospitime micro spring needle holder is a double spring instrument utilized for holding small needles in microsurgical procedures. It also features a locking mechanism.

PRODUCT TAGS

- Micro Spring Needle Holder Without Lock
- Ophthalmic Instruments
- Needle Holders
- Micro Surgery Instruments
- Needle Holder

TECHNICAL FOOTER

Product URL: <https://catalog.hospitime.com/product/micro-spring-needle-holder-without-lock-17-5-cm-straight/>
Document: Product Datasheet, Rev. 1 — Generated 2026-09-19
Manufacturer: Hospitime India — ISO 9001:2015 & CE Certified
Contact: info@hospitime.com | +91-11-4223-0205 | www.hospitime.com