




Live • Health • Life



Live • Health • Life

GENERAL SURGERY

ABOUT US

Background:

Hospitime is a specialised service provider in the Medico Surgical Field for more than six decades always intending to be a company that values Quality, Service and Commitment. Hospitime over its decades of operation has earned reputation as a quality manufacturer of Surgical Instruments and Tools. The company has been successfully servicing large number of Hospitals and Institutions in India and abroad with its manufactured products. Hospitime has also supplied surgical instruments and kits to various pilot programmes to NGO's, Govt. Institutions and Hospitals across India.

Dedication to Quality:

We employ the highest quality of raw materials. Hospitime is an ISO 9001:2008 certified company and all its products are CE compliance. This guarantee that our instruments are safe and fit for their intended use.


We owe our success to our prestigious customers who have maintained their trust in our products, over the years.

"Quality is our culture"

OUR CERTIFICATIONS



OUR MANUFACTURING PROCESS



WE SET THE BENCHMARK FOR PRODUCTION.

Did you know that ...?

- Metzenbaum scissors (as shown on Pg no. 05) are named after the American surgeon, Dr. Myron Firth Metzenbaum who designed the scissor to cut delicate tissues.

- Lister bandage scissor, (Pg no. 05) were invented by a British Surgeon Dr. Joseph Lister.

- Tungsten carbide or TC forceps and scissors (as shown on Pg no. 24) generally provide precise cut and long lasting grip than those without TC. They usually recognized by their gold-plated ring handles.

- Needle holder (Pg No. 09) and Hegar's dilator (Pg no.) both were developed by German gynecologist Alfred Hegar.

- Rongeurs are also called bone biters or bone nibblers, as their main function is to cut or remove small pieces of tissue or bone. Double action rongeurs provide much more power than single action, causing less hand fatigue.

- The primary use of a hemostat is to clamp and hold blood vessels so that the patient does not bleed to death.





No. 3
1001



No. 4
1004



No. 7
1007

Scalpel Handle, B.P



Straight
1012

Scissor, Fine (Pointed)
10.5cm



Straight
1018

Scissor, Kilner
11.5cm



Straight

1022	12.5cm
1024	15cm
1026	17.5cm
1028	20cm

1023	12.5cm
1025	15cm
1027	17.5cm
1029	20cm

Straight

Curved



Straight

1030	12.5cm
1032	15cm
1034	17.5cm
1036	20cm

1031	12.5cm
1033	15cm
1035	17.5cm
1037	20cm

Straight

Curved



Straight

1052	14.5cm
1054	16.5cm
1056	19cm
1058	21.5cm

1053	14.5cm
1055	16.5cm
1057	19cm
1059	21.5cm

Straight

Curved

Scissor, Dressing (Sharp X Sharp)

Scissor, Dressing (Sharp X Blunt)

Scissor, Mayo's



Scissor, Metzenbaum



Scissor, Suture Removal "HEATH"
15cm



Scissor, Gauge Cutting

1078	20cm
1079	12cm
1080	25cm



Scissor, Bandage Cutting "LISTER"



Scissor, Wire Cutting



Scissor, Pott's
17.5cm



Dissecting Forceps, Tissue



Dissecting Forceps, Gillies
15cm



Dissecting Forceps, Waugh
20cm



Dissecting Forceps, Dennis Browne
18cm




Dissecting Forceps, Russian



Dissecting Forceps, Jeweler
12cm



Dissecting Forceps, Suture Tying
10cm



Dissecting Forceps, Adsons
11cm



Hemostat Forceps, "MOSQUITO"



Hemostat Forceps, "SPENCERWELLS"



Hemostat Forceps, "KOCHER'S"



Hemostat Forceps, "MOYNIHAN"



Hemostat Forceps, Ligature "MIXTURE"



Hemostat Forceps, "LAHEY'S"
21cm



1272	10cm
1273	12.5cm
1274	15cm
1275	20cm
1276	25cm

Tissue Forceps, Allis



1283	10cm
1284	15cm
1285	20cm

Tissue Forceps, Babcock



1293	10cm
1294	12.5cm
1295	15cm

Towel Forceps, Backhaus



1296

Towel, Clip Cross Action
9cm



1301

Jar, Cheattle




1304

Sterilizer Forceps, Cheattle
25cm




1312	25cm
1313	30cm
1314	35cm
1315	40cm

Sterilizer Forceps, Bowl Lifting "HARRISON"




1321	15cm	Straight
1323	23cm	
1325	25cm	
1322	15cm	Curved
1324	23cm	
1326	25cm	

Intestinal Clamp, Doyen



1331	6cm	Blade
1332	10cm	
1333	15cm	

Intestinal Clamp, Payr's Crushing



Straight 1338

Intestinal Clamp, Lanes Twin



1351	10cm
1352	12.5cm
1353	15cm
1354	17.5cm
1355	20cm
1356	25cm

Needle Holder, Mayo Heagar



1360	15cm
1361	20cm

Needle Holder, Bozemann




1366

Needle Holder, Gillies



Straight 1377

Needle Holder, Micro Spring without Lock



Sharp
1392



Curved
1393

Retractor, Skin "KILNER"



2 Prongs
1397



3 Prongs
1398



4 Prongs
1399




5 Prongs
1400



6 Prongs
1401


Sharp



2 Prongs
1403



3 Prongs
1404



4 Prongs
1405



5 Prongs
1406



6 Prongs
1407

Blunt

Retractor, Volkman



1416

Retractor, Skin Hook



Sharp
1421



Blunt
1422

Retractor, Double Hook



1428	7cm
1429	12cm
1430	18cm
1431	25cm

Retractor, Langenbeck



1434

Retractor, Czerny



1438

Retractor, Allison Lung




1440

Retractor, Lane's Set of 2



1442

Retractor, Roux Set of 3



1447	25x1.2cm
1448	25x2.5cm
1449	30x1.2cm
1450	30x2.5cm
1451	30x3.8cm
1452	30x5cm
1453	30x6.3cm
1454	30x7.6cm

Length x Blade



1457	3.8cm
1458	5cm

Blade

Retractor, Kelly's



1460	5cm
1461	6.3cm
1462	7.6cm
1463	8.9cm
1464	10.1cm

Blade


Retractor, Morris



1468	5cm
1469	6.3cm
1470	7.6cm
1471	8.9cm
1472	10.1cm

Blade


Retractor, Doyen



1484

Self Retaining Retractor, Abdominal "BALFOUR"

Self Retaining Retractor, Rib's Spreader "FINOCHIETTO"



1495

Self Retaining Retractor, Pozzi with 3 blades

1496

Self Retaining Retractor, Millian Bladder with 5 blades

1497

Self Retaining Retractor, Alan-Parks



Self Retaining Retractor, Thyroid "JOLLS"



Self Retaining Retractor, Lateral Wall



Self Retaining Retractor, Gelpi

1509	15cm
1510	20cm



Self Retaining Retractor, Jansen
10cm



Self Retaining Retractor, Mollison



Self Retaining Retractor, Adson "BECKMANN"

1536	15cm
1537	20cm
1538	25cm
1539	30cm





1611

**Urethral Dilator, Lister
Set of 12**



1614

**Urethral Dilator, Clutton
Set of 12**



1619

Trocars, Suprapubic "NELSON"



1623

Trocars, Antrum "LICHTWITZ"



1628

Vein Strippers "VERICOSE"




1634

Clamp, Tubing



1636

Kidney Stone Forceps, Renal "RENDALLS"
Set of 4



3MM
HOSPITIME α

4MM
HOSPITIME α

5MM
HOSPITIME α

6MM
HOSPITIME α

7MM
HOSPITIME α

8MM
HOSPITIME α

9MM
HOSPITIME α

10MM
HOSPITIME α

11MM
HOSPITIME α

1639

Dilators, Gall Duct "BAKE" Set of 9




1642

Clamp, Satinsky



Straight
1646

Clamp, DeBakey Bull-Dog



Curved
1647



1652

Tissue Forceps, DeBakey Vascular





2168

Sthethoscope, Foetal Heart "PINARD"



2250	Small
2251	Medium
2252	Large

Vaginal Speculum, Kristeller



2258	Small
2259	Medium
2260	Large

Vaginal Speculum, Sims



2408	Small
2409	Medium
2410	Large

Vaginal Speculum, "CUSCO"



2416	Small
2417	Medium
2418	Large

Vaginal Speculum, "GRAVES"



2433

Retractor, Anterior Vaginal Wall



No. 1
2434



No. 2
2435



No. 3
2436



No. 4
2437



No. 5
2438



No. 6
2439

**Curette, Uterine Single Ended "SIMS"
Sharp/Blunt**



Sharp

2445	NO. 00
2446	NO. 0
2447	NO. 1
2448	NO. 2
2449	NO. 3
2450	NO. 4




Blunt

**Curette, Uterine Double Ended
(Sharp X Blunt)**



2477

**DILATORS, HEGAR
½MM
Set of 10**




2478

**DILATORS, HEGAR
Set of 8**



2512

Scissor, Episiotomy



2515

Scissor, Umbilical Cord (American Pattern)



2518

Speculum, Kogan Endocervical



2529

Biopsy Punch, Uterine



2534

Biopsy Punch, Cervical



2539

Biopsy Punch, Van Doren



2572	15cm
2573	20cm
2575	25cm

Forceps, Sponge Holding



2589

Uterine Forceps, Polypus "KELLY"



Curved
2593

Straight
2592

Uterine Forceps, Ovum



2628

Uterine Forceps, Vulsellum



2634

Uterine Forceps, Dressing "BOZEMANN"



2637

Uterine Manipulator



2641

Uterus Holding Forceps, "DARTIGUES"



2642

Uterus Holding Forceps, "SHIRODKAR'S"



2643

Uterus Holding Forceps, "SOMER"



2648

Forceps, Green Armytage



2652

Screw, Myoma



2658

Midwifery Forceps, Wrigley



Hysterectomy Forceps, Heaney Single Tooth



Hysterectomy Forceps, Heaney Double Teeth



Hysterectomy Forceps, "MAINGOT'S"



Hysterectomy Forceps, Rogers



Sound, Uterine "SIMS"




Female Catheter



Cannula, Hysto-Graphi with Lock Set of 3



Scissor, Metzenbaum



Dissecting Forceps, Tissue



Hemostat Forceps, "MOSQUITO"



Tissue Forceps, Allis



Tissue Forceps, Babcock




Needle Holder, Mayo Heagar



Scissor, Fine (Pointed)
10.5cm



Scissor, Kilner
11.5cm



Scissor, Mayo's



Scissor, Metzenbaum



Dissecting Forceps, Tissue



Dissecting Forceps, Adson
11.5cm



3573	15cm
3574	17.5cm
3575	20cm

Needle Holder, Mayo Heagar



3580	15cm
3581	17.5cm
3582	20cm

Needle Holder, Ryder



Curved
3597

Straight
3596

Needle Holder, Micro Spring without Lock



3608	15cm
3609	20cm

Tissue Forceps, Allis



3617	15cm
3619	20cm

Tissue Forceps, Babcock



3709

Cutter, Wire



Speculum, Heath Aural
Oval Shape



Speculum, Siegle Pneumatic
Set of 3



Tuning Forks, Hartmann

4065	128Hz
4066	256Hz
4067	512Hz
4068	1024Hz
4069	2048Hz
4070	4096Hz



Syringe, Ear



Knife, Sinus Sickle



Forceps, Hartmann Aural



Forceps, Micro Crocodile



4254

Suction Tube, Frazier



4267

Forceps, Tilley Nasal



4285

Scissor, Heymann Turbinectomy



Small
4230



Medium
4301



Large
4302

**Nasal Speculum, Thudichum
Set of 3**



4334

Nasal Speculum, Killian



4400

Forceps, Ash Septum



4524

Suction Cannula, Yankeur



4579

Depressor, Tongue
Set of 3



4623

Forceps, Tracheal Dilating "Trousseau"



4687

Curette, Adenoid




Small
4731

Medium
4732

Large
4733

Foceps, Magill Intratracheal Catheter



4835

Draffin Bipod Stand
With Plates



4900

**Mouth Gag, Boyle Davis
4 Tongue Plates**



4905

**Mouth Gag, Dingman
4 Tongue Plates**



4910	Child
4911	Adult

Mouth Gag, Doyan



4913

Mouth Gag, Heister



4915

Mouth Gag, Jennings



Straight
5138

Rongeur, Bone



Straight
5144

Forceps, Bone Cutting



5172

Self Centering Bone Holding Forceps



5267

Elevator, Periosteal



5312

Driver, Hexagonal Screw



5388

Hammer, Bone



5400

Bender, Plate



5446

Retractor, Hohmann



5521

Sleeve, Drill



5579

Guide, Drill



5630

Coupling Handle
T-Type




5684

Plaster Saw, Heavy Handle



5685

Plaster Saw, Engel




5745

Bone Drill, Universal



6027

**Gigle Saw, Handles
(Pair)**



6028

Gigle Saw, Wire



6083

Reny's Clip Applicator



6156

Hemostat Forceps, Dandy



6216


Punch Forceps, Disc



Upward
6216



Downward
6216



6256	1mm
6257	2mm
6258	3mm
6259	4mm
6260	5mm

Upward

Punch Forceps, Kerrison



Straight

Forceps, Bayonet

6289	17.5cm
6290	23cm



6347

Hook, Yasargil Spring Fish Type



6412

Brace, Hudson



6438

Retractor, Cloward With 10 Blades





7134

**Skin Grafting
Handle, Silver Skin**



7135

**Skin Grafting
Handle, Watson**



7136

**Skin Grafting
Handle, Humby**



908

**Finger Switch Hand Control
Electrosurgery Pencil**



Straight
922

Bipolar Forceps, Straight



Straight
923



Straight
939

Bipolar Forceps, Bayonet



Straight
940



Straight
955

**Bipolar Forceps, Bayonet
with Suction**




Straight
956



Straight
968

Cord, Bipolar Connection



Set of 4 Blades
980

Miller Laryngoscope



Set of 4 Blades
985

MacIntosh Conventional Laryngoscope



Lamp For Laryngoscope
999



Set of 4 Blades
990

MacIntosh Fiber Optic Laryngoscope



Set of 3 Blades
996

Flexi Tip Fiber Optic Laryngoscope



Breathing Bag

800	Paediatric
801	Child
802	Adult




Silicone Ventouse Cup

805	Small
806	Medium
807	Large



Cold Light Source
812




Fibre Optic Cable

814	230x3.5mm	LxT
815	230x4.5mm	



Fibre Optic Head Light Band
818



Formalin Chamber

825	14"x7"x6"	LxBxH
826	20"x8"x8"	
827	23"x8"x8"	
828	26"x8"x8"	



**Sterilization Box
with Silicone Mat**

832	8"x5"x1"	LxBxH
833	9"x8"x1"	



Instrument Tray

855	8"X6"	860	12"X8"
856	9"X5"	861	12"X10"
857	9"X6"	862	14"X10"
858	10"X8"	863	15"X12"
859	11"X7"	864	18"X12"



Dressing Drum

870	6"X6"
871	9"X9"
872	11"X9"
873	14"X9.5"
874	15"X12"



Autoclave
14"X12"
878



Bed Pen

888	Male
889	Female



Basin

894	30cm
895	35cm
896	40cm




Kidney Tray

897	15cm
898	20cm
899	25cm
900	30cm




Bone Drill Machine
5747



Plaster Saw Cutter Electric
5749



INSTRUMENT SET

S.No	Surgery	Name	Product Code
1.	General	Basic Laparotomy Set	650
2.	General	Minor Basic Instrument Set	653
3.	General	Plastic Set	658
4.	General	Hand Set	659
5.	General	Pediatric Laparotomy Set	660
6.	General	Appendectomy & Hernia Set	661
7.	General	Common Duct Instruments	662
8.	General	Bowel Resection Set	663
9.	General	Abdominoperineal Resection Set	664
10.	General	Thyroidectomy Set	665
11.	Genito-Urinary	Nephrectomy Set	670
12.	Genito-Urinary	Pyeloplasty & Ureteroplasty Set	671
13.	Genito-Urinary	Suprapubic & Retropubic Prostatectomy Set	672
14.	Genito-Urinary	Vasectomy-Meatotomy Set	673
15.	Genito-Urinary	Tuboplasty Set	674
16.	Rectal	Rectal Instrument Set	680
17.	Rectal	Pilonidal Cyst Excision Set	681
18.	Rectal	Hemorrhoidectomy Set	682
19.	Ob/Gyn	Dilatation & Curettage Set	686
20.	Ob/Gyn	Basic Vaginal Instrument Set	687
21.	Ob/Gyn	Abdominal Gynecological Set	688
22.	Ob/Gyn	Laparoscopy Set	689
23.	Ob/Gyn	Bartholin Cyst Excision Set	690
24.	Ob/Gyn	Vaginal Tubal Ligation Set	691
25.	ENT	Basic Ear Set	697
26.	ENT	Myringotomy Set	698
27.	ENT	Basic Nasal Set	699
28.	ENT	Minor Nasal Set	700
29.	ENT	Major Nasal Set	701
30.	ENT	Tympanoplasty	702
31.	ENT	Mastoid Surgery Set	703
32.	ENT	Septorhinoplasty Set	704
33.	ENT	Paranasal Sinus Surgery Set	705
34.	ENT	DCR Surgery Set	706
35.	ENT	FESS Surgery Set	707
36.	ENT	Tonsillectomy & Adenoidectomy Set	708
37.	ENT	Tracheostomy Set	709
38.	Orthopedic	Basic Major Orthopedic Set	715
39.	Orthopedic	Basic Orthopedic Soft Tissue Set	716
40.	Orthopedic	Basic Total Hip Set	717
41.	Orthopedic	Knee Set	718
42.	Orthopedic	Amputation Set	719
43.	Neuro/Spine	Basic Craniotomy Set	726
44.	Neuro/Spine	Ventriculography Set	727
45.	Neuro/Spine	Neurosurgical - Laminectomy Set	728
46.	Neuro/Spine	Skull Traction Set	729
47.	Thoracic & Cardiac	Basic Instrument Set	734
48.	Thoracic & Cardiac	Thoracotomy Set	735
49.	Thoracic & Cardiac	Closed Thoracostomy Set	736
50.	Thoracic & Cardiac	Open Thoracostomy Set	737
51.	Thoracic & Cardiac	Lobectomy & Segmental Lung Resection Set	738
52.	Thoracic & Cardiac	Basic Open Heart Set	739
53.	Thoracic & Cardiac	Major Cardiovascular Set	740
54.	Thoracic & Cardiac	Aortic-Corony Graft Set	741
55.	Ophthalmic	Basic Eye Set	746
56.	Ophthalmic	Cataract Surgery Set	747
57.	Ophthalmic	Eye Muscle & Enucleation Set	748
58.	Ophthalmic	Retinal Surgery Set	749
59.	Ophthalmic	Dacryocystorhinostomy Set	750
60.	Suture	E.R. Suture Set	757
61.	Suture	Small Suture Set	758



Alphabetical Index

ADENOID CURETTE.....	25	FLEXI TIP LARYNGOSCOPE.....	32	MOLLISON SELF RETAINING RETRACTOR.....	10
ADSONS DISSECTING FORCEPS.....	03	FOETAL HEART "PINARD" STHETHOSCOPE.....	14	MORRIS RETRACOR.....	08
ADSONS DISSECTING FORCEPS TUNGSTEN CARBIDE.....	21	FORMALIN CHAMBER.....	33	MOSQUITO HEMOSTAT FORCEPS.....	04
ALAN-PARKS SELF RETAINING RETRACTOR.....	09	FRAZIER SUCTION TUBE.....	24	MOSQUITO HEMOSTAT FORCEPS TITANIUM.....	20
ALLIS TISSUE FORCEPS.....	05	GALL DUCT "BAKES" DILATOR.....	13	MOYNIHAN HEMOSTAT FORCEPS.....	04
ALLIS TISSUE FORCEPS TITANIUM.....	20	GAUGE CUTTING SCISSOR.....	02	OVUM FORCEPS.....	17
ALLIS TISSUE FORCEPS TUNGSTEN CARBIDE.....	22	GELPI SELF RETAINING RETRACTOR.....	10	PAYR'S CRUSHING INTESTINAL CLAMP.....	06
ALLISON LUNG RETRACTOR.....	08	GIGLE SAW HANDLES.....	29	PERIOSTEAL ELEVATOR.....	27
AMERICAN PATTERN UMBILICAL SCISSOR.....	16	GIGLE SAW WIRE.....	29	PILE BAND RUBBER.....	11
ANTERIOR VAGINAL WALL RETRACTOR.....	15	GILLIES DISSECTING FORCEPS.....	03	PLASTER SAW CUTTER ELECTRIC.....	34
ANTRUM TROCAR LICHTWITZ.....	12	GILLIES NEEDLE HOLDER WITH SCISSOR.....	06	PLATE BENDER.....	27
ASH SEPTUM FORCEPS.....	24	GRAVES VAGINAL SPECULUM.....	14	POLYPUS "KELLY" UTERINE FORCEPS.....	17
AUTOCLAVE.....	34	GREEN ARMYTAGE FORCEPS.....	18	POTT'S SCISSOR.....	02
BABCOCK TISSUE FORCEPS.....	05	HARRISON BOWL LIFTING STERILIZER FORCEPS.....	05	REN'S CLIP APPLICATOR.....	29
BABCOCK TISSUE FORCEPS TITANIUM.....	20	HARTMANN AURAL FORCEPS.....	23	RETRACTOR COLLIN'S (POZZI) SELF RETAINING.....	09
BABCOCK TISSUE FORCEPS TUNGSTEN CARBIDE.....	22	HARTMANN TUNING FORK.....	23	ROGERS HYSTERECTOMY CLAMP.....	19
BALFOUR SELF RETAINING RETRACTOR.....	09	HEANEY'S DOUBLE TEETH HYSTERECTOMY CLAMP.....	19	ROUX RETRACTOR.....	08
BANDAGE SCISSOR "LISTER".....	02	HEANEY'S SINGLE TOOTH HYSTERECTOMY CLAMP.....	19	RUSSIAN DISSECTING FORCEPS.....	03
BASIN.....	34	HEATH AURAL SPECULUM.....	23	RYDER NEEDLE HOLDER TUNGSTEN CARBIDE.....	22
BAYONET BIPOLAR FORCEPS.....	31	HEAVY HANDLE PLASTER SAW.....	28	SATINSKY CLAMP.....	13
BAYONET FORCEPS.....	30	HEGARS UTERINE DILATOR DOUBLE ENDED.....	15	SCALPLE HANDLE B.P.....	01
BAYONET WITH SUCTION BIPOLAR FORCEPS.....	31	HEISTER MOUTH GAG.....	26	SELF CENTERING BONE HOLDING FORCEPS.....	27
BECKMANN SELF RETAINING RETRACTOR.....	10	HEMORRHOID FORCEPS.....	11	SHIRODKAR'S UTERUS HOLDING FORCEPS.....	18
BED PEN FEMALE.....	34	HEMORRHOID LIGATING FORCEPS.....	11	SIEGLE PNEUMATIC SPECULUM.....	23
BIPOLAR CONNECTION CORD.....	31	HEXAGONAL SCREW DRIVER.....	27	SILICONE VENTEUSE CUP.....	33
BONE CUTTING FORCEPS.....	27	HEYMANN TURBINECTOMY SCISSOR.....	24	SILVER SKIN GRAFTING HANDLE.....	31
BONE DRILL MACHINE.....	34	HOHMANN RETRACTOR.....	28	SIMS RECTAL SPECULUM.....	11
BONE HAMMER.....	27	HUDSON BRACE.....	30	SIMS UTERINE SOUND.....	19
BONE ROUNGER.....	27	HUMBY SKIN GRAFTING HANDLE.....	31	SIMS VAGINAL SPECULUM.....	14
BOYLE DAVIS MOUTH GAG.....	26	HYSTO-GRAPHI CANNULA WITH LOCK.....	19	SINUS SICKLE KNIFE.....	23
BOZEMANN NEEDLE HOLDER.....	06	INSTRUMENT TRAY.....	34	SKIN HOOK RETRACTOR.....	07
BREATHING BAG.....	33	JANSEN SELF RETAINING RETRACTOR.....	10	SOMER'S UTERUS HOLDING FORCEPS.....	18
BULL DOG CLAMP ATRAUMATIC.....	13	JENNINGS MOUTH GAG.....	26	SPENCERWELLS HEMOSTAT FORCEPS.....	04
CERVICAL BIOPSY PUNCH FORCEPS.....	16	JEWELER DISSECTING FORCEPS.....	03	SPONGE HOLDING FORCEPS CURVED.....	17
CHEATLE JAR.....	05	KELLY RECTAL SPECULUM.....	11	STERILIZATION BOX WITH SILICONE MAT.....	33
CHEATLE STERILIZER FORCEPS.....	05	KELLY'S RETRACTOR.....	08	STRAIGHT BIPOLAR FORCEPS.....	31
CLOWARD RETRACTOR.....	30	KERRISON PUNCH FORCEPS.....	30	SUTURE CUTTING SCISSOR "HEATH".....	02
COLD LIGHT SOURCE.....	33	KIDNEY STONE FORCEPS RANDALL.....	13	SUTURE TYING DISSECTING FORCEPS.....	03
COUPLING HANDLE T-TYPE.....	28	KIDNEY TRAY.....	34	THUDICHUM.....	24
CUSCO VAGINAL SPECULUM.....	14	KILLIAN SPECULUM.....	24	THYROID "JOLLS" SELF RETAINING RETRACTOR.....	10
CZERNY RETRACTOR.....	07	KILNER SCISSOR.....	01	TILLEY NASAL FORCEPS.....	24
DANDY HEMOSTAT FORCEPS.....	29	KILNER SCISSOR TUNGSTEN CARBIDE.....	21	TISSUE DISSECTING FORCEPS.....	03
DARTIGUES UTERUS HOLDING FORCEPS.....	18	KILNER SKIN RETRACTOR.....	07	TISSUE DISSECTING FORCEPS TITANIUM.....	20
DEAVERS RETRACTOR.....	08	KOCHER'S HEMOSTAT FORCEPS.....	04	TISSUE DISSECTING FORCEPS TUNGSTEN CARBIDE.....	21
DeBakey VASCULAR TISSUE FORCEPS.....	13	KOGAN ENDOCERVICAL SPECULUM.....	16	TONGUE DEPRESSOR.....	25
DENNIS BROWNE DISSECTING FORCEPS.....	03	KRISTELLER VAGINAL SPECULUM.....	14	TONSIL/PILE NEEDLE.....	11
DINGMAN MOUTH GAG.....	26	LAHEY'S HEMOSTAT FORCEPS.....	04	TOWEL CLIP CROSS ACTION.....	05
DISC PUNCH FORCEPS.....	29	LANE'S RETRACTOR.....	08	TOWEL FORCEPS "BACKHAUS".....	05
DOUBLE HOOK RETRACTOR SHARP & BLUNT.....	07	LANE'S TWIN INTESTINAL CLAMP.....	06	TRACHEAL DILATING "TROUSSEAU" FORCEPS.....	25
DOYAN MOUTH GAG.....	26	LANGENBECK RETRACTOR.....	07	TROCAR & CANULA.....	12
DOYEN INTESTINAL CLAMP.....	06	LATERAL WALL SELF RETAINING RETRACTOR.....	10	TUBING CLAMP.....	12
DOYEN MYOMA SCREW.....	18	LIGATURE HEMOSTAT FORCEPS "MIXTER".....	04	UNIVERSAL BONE DRILL.....	28
DOYEN RETRACTOR.....	08	MAGILL INTRATRACHEAL CATHETER FOCEPS.....	25	URETHRAL BOUGIES CLUTTON.....	12
DRESSING "BOZEMANN" UTERINE FORCEPS.....	17	MAINGOT HYSTERECTOMY CLAMP.....	19	URETHRAL BOUGIES LISTER.....	12
DRESSING DRUM.....	34	MAYO HEGAR NEEDLE HOLDER.....	06	UTERINE BIOPSY PUNCH FORCEPS.....	16
DRESSING SCISSOR.....	01	MAYO HEGAR NEEDLE HOLDER TITANIUM.....	20	UTERINE CURETTE DOUBLE ENDED.....	15
DRILL GUIDE.....	28	MAYO HEGAR NEEDLE HOLDER TUNGSTEN CARBIDE.....	22	UTERINE CURETTE SINGLE SHARP.....	15
DRILL SLEEVE.....	28	MAYO'S SCISSOR.....	01	UTERINE MANIPULATOR.....	17
EAR SYRINGE.....	23	MAYO'S SCISSOR TUNGSTEN CARBIDE.....	21	VAN DOREN BIOPSY PUNCH FORCEPS.....	16
EARDRAFFIN BIPOD STAND.....	25	McIntosh CONVENTIONAL LARYNGOSCOPE.....	32	VIEN STRIPPER VERICOSE.....	12
ENGEL PLASTER SAW.....	28	McIntosh FIBER OPTIC LARYNGOSCOPE.....	32	VOLKMAN RETRACTOR 2 PRONG SHARP & BLUNT.....	07
EPISIOTOMY SCISSOR.....	16	METZENBAUM SCISSOR.....	02	VULSELLUM FORCEPS.....	17
FIBRE OPTIC CABLE.....	33	METZENBAUM SCISSOR TITANIUM.....	20	WATSON SKIN GRAFTING HANDLE.....	31
FIBRE OPTIC HEAD LIGHT BAND.....	33	METZENBAUM SCISSOR TUNGSTEN CARBIDE.....	21	WAUGH DISSECTING FORCEPS.....	03
FINE SCISSOR POINTED SHARP.....	01	MICRO CROCODILE FORCEPS.....	23	WIRE CUTTER 7" TUNGSTEN CARBIDE.....	22
FINE SCISSOR POINTED SHARP TUNGSTEN CARBIDE.....	21	MICRO SPRING NEEDLE HOLDER TUNGSTEN CARBIDE.....	22	WIRE CUTTING SCISSOR.....	02
FINGER SWITCH HAND CONTROL.....	31	MICRO SPRING NEEDLE HOLDER WITH LOCK.....	06	WRIGLEY MIDWIFERY FORCEPS.....	18
FINOCHIETTO SELF RETAINING RETRACTOR.....	09	MILLER LARYNGOSCOPE.....	32	YANKEUR SUCTION CANNULA.....	25
FISTULA PROBE.....	11	MILLIN BLADDER SELF RETAINING RETRACTOR.....	09	YASARGIL SPRING FISH TYPE.....	30



Clean and lubricate the brand new Hospitime instruments before the first sterilization. After use, distinct, clean, lubricate and sterilize the instruments as soon as possible.

Disinfect chemothermally, thermally or by immersing into a combined disinfecting & cleaning solution. Delicate instruments can be cleaned by ultrasonic treatment.

Rinse thoroughly and carefully after cleaning if possible, use demineralised water. Dry sufficiently after rinsing

Treat the instruments with joints, ratchets or threads with lubricating agent based on paraffin oil. We strongly recommend to follow the sterilization instructions of the manufacturer.

Prior to sterilization, only close the first ratchet on the instruments. Normally, autoclaving is performed with saturated water stream at 134°C.

When instruments are hot-air sterilized, the temperature should be between 180°C and 200°C. Corroded instruments should be discarded immediately



OB/GYN Instruments



Neuro/Spine Instruments



Orthopedic Instruments



ENT Instruments



Ophthalmic Instruments



Plastic Surgery Instruments



Diagnostic Products



Tungsten Carbide Instruments



Titanium Instruments



Hollowares and Utensils



HOSPITIME INDIA

Manufacturer & Exporter of Surgical Instruments

Office: BM-32 East, Shalimar Bagh, DELHI-110088 (DL)

Works: 1538, HSIIDC, Ind. Estate, Rai, Dist. SONPAT (HR)

Telephone: +91-11-27498687, 42230205

Email: info@hospitime.com | www.hospitime.com

