
ADSON DISSECTING FORCEPS

11cm 1X2 Teeth | SKU 1149



PRODUCT OVERVIEW

Summary	1 x 2 teeth. Overall length 4-1/4in (…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Adsons Instrument Type: Tissue Forceps Curvature: Straight Length: 11 cm Tip Type: 1x2 Teeth (Tooth) Handle: Thumb Material: Stainless Steel Alias Names: adsons tissue forceps, adson dressing forceps Surgical Instrument: Hospitime adsons dissecting forceps 11cm (tooth) are for holding and controlling delicate tissues. This forceps have a wide, flat thumb region that is commonly serrated. The jaws are short and the tips are narrow.

PRODUCT TAGS

Adson Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Plastic / Dermatology Instruments | Forceps & Clamps

TECHNICAL FOOTER

Product URL: https://catalog.hospitime.com/product/adson-dissecting-forceps-11-cm-1x2-teeth/
Document: Product Datasheet, Rev. 1 — Generated 2026-11-03
Manufacturer: Hospitime India — ISO 9001:2015 & CE Certified
Contact: info@hospitime.com | +91-11-4223-0205 | www.hospitime.com