
ADSON DISSECTING FORCEPS

11cm Plain | SKU 1148



PRODUCT OVERVIEW

Summary	Serrated. Overall length 4-1/4in (11…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Adsons Instrument Type: Tissue Forceps Curvature: Straight Length: 11 cm Tip Type: Plain Handle: Thumb Material: Stainless Steel Alias Names: adsons tissue forceps, adson dressing forceps Surgical Instrument: Hospitime adsons dissecting forceps 11cm (serrated) are for holding and controlling delicate tissues. This forceps have a wide, flat thumb region that is commonly serrated. The jaws are short and the tips are narrow.

PRODUCT TAGS

Adson Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Plastic / Dermatology Instruments | Forceps & Clamps