
ALLIS TISSUE FORCEPS

17.5cm | SKU 1275



PRODUCT OVERVIEW

Summary	4 x 5 teeth. Overall length 7in (17.5…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Tissue holding forceps Length: 17.5 cm Tip Type: 4x5 teeth Handle: Ratched, Finger Rings Material: Stainless Steel Alias Names: Allice forceps, Tissue grasping forceps Surgical Instrument: Hopsitime allis tissue forceps are ratcheted forceps, utilized for grasping organs and elusive or thick tissue surgery. They are regularly utilized with breast, tonsil, thyroid and vaginal tissues.

PRODUCT TAGS

Allis Tissue Forceps | General Surgery Instruments | Forceps, Clamps & Clips