BABCOCK TISSUE FORCEPS

20cm | SKU 1285

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 8in (20cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Tissue holding forceps Length: 20 cm Tip Type: Circumferential with serrated jaws (rounded grasping surface) Handle: Ratched, Finger Rings Material: Stainless Steel Alias Names: Tissue grasping forceps Surgical Instrument: Hospitime Babcock tissue forceps are ratcheted forceps utilized to handle sensitive tissues. They are often utilized with intestinal and laparotomy techniques. Babcock Forceps are similar to Allis forceps. The jaws are circumferential, and the tips are serrated.

PRODUCT TAGS

Babcock Tissue Forceps | General Surgery Instruments | Forceps, Clamps & Clips