BIRKETT HEMOSTAT FORCEPS

19cm Straight | SKU 1260

[product image not available]

PRODUCT OVERVIEW

Summary	Straight. Overall length 7-1/2in (19…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Birkett Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 19 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Birkett artery forceps Surgical Instrument: Hospitime Birkett Hemostatic Forceps are slim, ratcheted, finger ring forceps utilized for clamping vessels and controlling blood float. they're similar to the Moynihan hemostat, with completely serrated jaws and commonly utilized in surgical procedures inside the abdominal area; nonetheless, the Birkett hemostat has somewhat bigger jaws and is marginally longer. This instrument is straight with an overall length of 7Â½ inches.

PRODUCT TAGS

Birkett Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips