BOZEMANN NEEDLE HOLDER

20cm | SKU 1361

[product image not available]

PRODUCT OVERVIEW

Summary	Bayonet. Overall length 8in (20cm).
Category	General Surgery Instruments, Needl…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Bozeman-Heaney Length: 20 cm Instrument Type: Needle Holder Curvature: 'S' Curved Blade Type: Cross-serrated with groove Tip Type: Blunt Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Bozeman-Heaney needle holder, Bozeman needle holder, Curved needle holder, Bayonet needle holder Surgical Instrument: Hospitime Bozeman needle holder is a ratcheted instrument, used in multiple types of surgical techniques to hold a suturing needle for closing wounds during suturing and surgical procedures. The shanks of this instrument are slightly curved and therefore provide greater visibility.

PRODUCT TAGS

Bozemann Needle Holder | General Surgery Instruments | Needle Holder