CHEATLE STERILIZER FORCEPS HEAVY PATTERN

28cm | SKU 1305

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 11in (28cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: LHP Instrument Type: Sterilization forceps Length: 28 cm Tip Type: Plain (No Serration) Handle: Non-Ratcheted Material: Stainless Steel Surgical Instrument: Hospitime cheatle sterilization forceps utilized for placing and removing instruments within a disinfection container.

PRODUCT TAGS

Cheatle Sterilizer Forceps Heavy Pattern (LHP) | General Surgery Instruments | Forceps, Clamps & Clips