CHEATLE STERILIZER FORCEPS WITHOUT SERRA

25cm | SKU 1304

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 10in (25cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Sterilization forceps Length: 25 cm Tip Type: Plain (No Serration) Handle: Non-Ratcheted & Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime cheatle sterilization forceps is a non-ratcheted finger ring forceps, utilized for placing and removing instruments within a disinfection container.

PRODUCT TAGS

Cheatle Sterilizer Forceps without Serration | General Surgery Instruments | Forceps, Clamps & Clips