
CRILE HEMOSTAT FORCEPS

15cm Straight | SKU 1188



PRODUCT OVERVIEW

Summary	Straight. Overall length 6in (15cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Crile Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 15 cm Tip Type: Fine & Horizontally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Crile artery forceps Surgical Instrument: Hospitime crile hemostat forceps is a ratcheted forceps utilized for clipping veins or tissue before cauterization or ligation. Jaws are half the length of the shank make this instrument remarkable.

PRODUCT TAGS

Crile Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips