CZERNY RETRACTOR

20cm | SKU 1435

[product image not available]

PRODUCT OVERVIEW

Summary	Double ended. Overall legth 8in (20…
Category	General Surgery Instruments, Retra…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Czerny Instrument Type: Retractors Blade Type: Right Angled- Flat blade & 2 spikes Length: 20 cm Tip Type: Curved Ends Handle: Flat handle with cut Material: Stainless Steel Alias Names: Czerny rectal retractor Surgical Instrument: Hospitime czerny retractor has two spikes toward one side and a broad surface on the opposite. It is commonly utilised during closure of a midline laparotomy to retract the higher or lower edges than supporting the inclusion of the primary mass closure stitch.

PRODUCT TAGS

Czerny Retractor | General Surgery Instruments | Retractor