DENNIS BROWNE DISSECTING FORCEPS

18cm | SKU 1120



PRODUCT OVERVIEW

Summary	Flat handle. Ring tips. Overall length…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Dennis-Browne Instrument Type: Tissue Forceps Curvature: Straight Length: 18 cm Tip Type: Cup-shaped with holes Handle: Thumb Material: Stainless Steel Surgical Instrument: Hospitime dennis browne dissecting forceps mostly used to hold the tonsil during tonsillectomy procedures, facilitate dissection and avoid harm to structures in the bed of the tonsil.

PRODUCT TAGS

Dennis Browne Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | ENT Instruments | Tonsillectomy