DOUBLE HOOK RETRACTOR

SKU 1421

[product image not available]

PRODUCT OVERVIEW

Summary	Double hook. Sharp or Blunt.
Category	General Surgery Instruments, Retra…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Retractors Prongs: 2 prong Tip Type: Sharp or Blunt Handle: Round Handle Material: Stainless Steel Surgical Instrument: Hospitime Double Hook contains two sharp prongs toward the end of a straight rod that contains a flat handle for enhanced grip. The skin hook is generally utilized s a part of surgeries where sensitive skin needs to be retracted, for example, retracting skin of the trachea. Depending upon the sort of skin or the kind of surgical procedure.

PRODUCT TAGS

Double Hook Retractor | General Surgery Instruments | Retractor