DOYEN INTESTINAL CLAMP

15cm Straight | SKU 1321

[product image not available]

PRODUCT OVERVIEW

Summary	Straight. Overall length 6in (15cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Doyen Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 15 cm Tip Type: Longitudinally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Intestinal occlusion clamp, Doyen forceps Surgical Instrument: Hospitime doyen intestinal clamp is a ratcheted & non-crushing clamp utilized to occlude the lumen of the intestine and the blood vessels as they are specially designed to handle delicate structures, large organs, for example, lungs or digestion systems. The Doyen forceps have long jaws and shorter shanks.

PRODUCT TAGS

Doyen Intestinal Clamp | General Surgery Instruments | Forceps, Clamps & Clips