DOYEN INTESTINAL CLAMP

23cm Curved | SKU 1324

[product image not available]

PRODUCT OVERVIEW

Summary	Curved. Overall length 9in (23cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Doyen Instrument Type: Hemostatic Clamps Blade Type: Curved Length: 23 cm Tip Type: Longitudinally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Kocher Intestinal Clamp Atraumatic, Intestinal occlusion clamp Surgical Instrument: Hospitime doyen intestinal clamp is a ratcheted forceps utilized to occlude the lumen of the intestine and the blood vessels as they are specially designed to handle delicate structures.

PRODUCT TAGS

Doyen Intestinal Clamp | General Surgery Instruments | Forceps, Clamps & Clips