# DOYEN INTESTINAL CLAMP

25cm Curved | SKU 1326

[ product image not available ]

## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Curved. Overall length 10in (25cm). |
| **Category** | General Surgery Instruments, Force… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Doyen Instrument Type: Hemostatic Clamps  Blade Type: Curved Length: 25 cm Tip Type: Longitudinally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Kocher Intestinal Clamp Atraumatic, Intestinal occlusion clamp Surgical Instrument: Hospitime doyen intestinal clamp is a ratcheted forceps utilized to occlude the lumen of the intestine and the blood vessels as they are specially designed to handle delicate structures.

## PRODUCT TAGS

Doyen Intestinal Clamp | General Surgery Instruments | Forceps, Clamps & Clips