DOYEN INTESTINAL CLAMP

25cm Straight | SKU 1325

[product image not available]

PRODUCT OVERVIEW

Summary	Straight. Overall length 10in (25cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Doyen Instrument Type: Hemostatic Clamps Blade Type: Straight Length: 15 cm Tip Type: Longitudinally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Kocher Intestinal Clamp Atraumatic, Intestinal occlusion clamp Surgical Instrument: Hospitime doyen intestinal clamp is a ratcheted forceps utilized to occlude the lumen of the intestine and the blood vessels as they are specially designed to handle delicate structures.

PRODUCT TAGS

Doyen Intestinal Clamp | General Surgery Instruments | Forceps, Clamps & Clips