# DRESSING SCISSORS (BLUNT X BLUNT)

15cm Curved | SKU 1041




## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Blunt/Blunt. Curved. Overall length 6… |
| **Category** | General Surgery Instruments, Scisso… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Dressing Instrument Type: Scissors Blade Type: Curved Length: 15 cm Tip Type: Blunt/Blunt Handle: Finger Rings Material: Stainless Steel Alias Names: Operating scissors, Surgical scissors Surgical Instrument: Hospitime dressing scissors have three accessible tip designs (sharp/sharp, sharp/blunt & blunt/blunt). These scissors are used for general, multipurpose cutting and dissecting. Sharp dressing scissors cut tougher tissue and are as often as possible utilized for dissecting, whereas blunt dressing scissors are as often as possible utilized cutting materials such as sutures, tubing and dressings or gauzes.

## PRODUCT TAGS

Dressing Scissors (Blunt x Blunt) | General Surgery Instruments | Scissors