DRESSING SCISSORS (BLUNT X BLUNT)

17.5cm Straight | SKU 1042



PRODUCT OVERVIEW

Summary	Blunt/Blunt. Straight. Overall length…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Dressing Instrument Type: Scissors Blade Type: Straight Length: 17.5 cm Tip Type: Blunt/Blunt Handle: Finger Rings Material: Stainless Steel Alias Names: Operating scissors, Surgical scissors Surgical Instrument: Hospitime dressing scissors have three accessible tip designs (sharp/sharp, sharp/blunt & blunt/blunt). These scissors are used for general, multipurpose cutting and dissecting. Sharp dressing scissors cut tougher tissue and are as often as possible utilized for dissecting, whereas blunt dressing scissors are as often as possible utilized cutting materials such as sutures, tubing and dressings or gauzes.

PRODUCT TAGS

Dressing Scissors (Blunt x Blunt) | General Surgery Instruments | Scissors