DRESSING SCISSORS (SHARP X SHARP)

17.5cm Curved | SKU 1027



PRODUCT OVERVIEW

Summary	Sharp/Sharp. Curved. Overall lengt…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Dressing Instrument Type: Scissors Blade Type: Curved Length: 17.5 cm Tip Type: Sharp/Sharp Handle: Finger Rings Material: Stainless Steel Alias Names: Operating scissors, Surgical scissors Surgical Instrument: Hospitime dressing scissors have three accessible tip designs (sharp/sharp, sharp/blunt & blunt/blunt). These scissors are used for general, multipurpose cutting and dissecting. Sharp dressing scissors cut tougher tissue and are as often as possible utilized for dissecting, whereas blunt dressing scissors are as often as possible utilized cutting materials such as sutures, tubing and dressings or gauzes.

PRODUCT TAGS

Dressing Scissors (Sharp x Sharp) | General Surgery Instruments | Scissors