GILLIES DISSECTING FORCEPS

15cm 1X2 Teeth | SKU 1116



PRODUCT OVERVIEW

Summary	Fine jaws and 1 x 2 teeth. Overall le…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Gillies Instrument Type: Tissue Forceps Curvature: Straight Length: 15 cm Tip Type: Fine tip with 1x2 Teeth (Tooth) Handle: Thumb Material: Stainless Steel Alias Names: Thumb forceps, kocher forceps, tooth forceps, fine tissue forceps Surgical Instrument: Hospitime gillies forceps are thumb forceps frequently utilized in small surgical locales to handle tissue. They have fine, contract jaws that are normally utilized in plastic surgery procedures.

PRODUCT TAGS

Gillies Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Plastic / Dermatology Instruments | Forceps & Clamps