GILLIES DISSECTING FORCEPS

15cm Plain | SKU 1115



PRODUCT OVERVIEW

Summary	Fine and serrated jaws. Overall leng…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Gillies Instrument Type: Tissue Forceps Curvature: Straight Length: 15 cm Tip Type: Fine tip with serrated jaws (Plain) Handle: Thumb Material: Stainless Steel Alias Names: Thumb forceps, Fine tissue forceps Surgical Instrument: Hospitime gillies forceps are thumb forceps frequently utilized in small surgical locales to handle tissue. They have fine, contract jaws that are normally utilized in plastic surgery procedures.

PRODUCT TAGS

Gillies Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Plastic / Dermatology Instruments | Forceps & Clamps