GILLIES NEEDELE HOLDER WITH SCISSORS

15cm | SKU 1366

[product image not available]

PRODUCT OVERVIEW

Summary	Combination needle holder and scis…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Gillies Length: 15 cm Instrument Type: Needle Holder Curvature: Curved Blade Type: Cross-serrated with groove & sharp scissor blade Tip Type: Blunt Handle: Finger Ring/offset ring Material: Stainless Steel Alias Names: Sheehan-Gillies needle holder Surgical Instrument: Hospitime gillies needle holder is a is a combination of the needle holder and scissors, used in multiple types of surgical techniques to hold a suturing needle for closing wounds during suturing and surgical procedures.

PRODUCT TAGS

Gillies Needele Holder With Scissors | General Surgery Instruments | Scissors | Needle Holder