IRIS SCISSORS FINE,

10.5cm Conjuctival Curved | SKU 1015



PRODUCT OVERVIEW

Summary	Sharp/Sharp. Curved. Overall length…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Iris Instrument Type: Scissors Blade Type: Curved Length: 10.5 cm (4¼") Handle: Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime iris scissors (curved) is a cutting instrument and are very small, with an extremely sharp tip (tiny like the iris in your eye.).

PRODUCT TAGS

Iris Scissors Fine, | General Surgery Instruments | Scissors