JEWELER DISSECTING FORCEPS

12cm Curved | SKU 1127



PRODUCT OVERVIEW

Summary	Curved. Overall length 4-3/4in (12cm…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Jeweler Instrument Type: Tissue Forceps Curvature: Curved Length: 12 cm Tip Type: Pointed & Fine Handle: Flat Alias Names: Micro forceps (jeweler type) Material: Stainless Steel Surgical Instrument: Hospitime jeweler forceps 12cm (curved) are grasping forceps utilized in microsurgical procedures. They have flat handles and fine & pointed tips to handle fine tissues and materials.

PRODUCT TAGS

Jeweler Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips