KILNER NEEDLE HOLDER

13cm | SKU 1362

[product image not available]

PRODUCT OVERVIEW

Summary	Bayonet. Overall length 5-1/4in (13…
Category	General Surgery Instruments, Needl…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Kilner Length: 13 cm Instrument Type: Needle Holder Curvature: 'S' Curved Blade Type: Cross-serrated with groove Tip Type: Blunt Handle: Ratchet & Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime Kilner needle holder is a ratcheted instrument, used in multiple types of surgical techniques designed to hold a fine suturing needle for closing wounds during suturing and surgical procedures. Features cross-hatched jaws to maintain the grip on needles.

PRODUCT TAGS

Kilner Needle Holder | General Surgery Instruments | Needle Holder