KILNER SCISSORS

11.5cm Straight | SKU 1018



PRODUCT OVERVIEW

Summary	Blunt/Blunt. Straight. Overall length…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Kilner Instrument Type: Scissors Blade Type: Straight Length: 11.5 cm Handle: Finger Rings Material: Stainless Steel Alias Names: Knapp scissor Surgical Instrument: Hospitime kilner scissors (straight) designed for cutting delicate tissue and blunt dissection.

PRODUCT TAGS

Kilner Scissors | General Surgery Instruments | Scissors