KOCHER'S HEMOSTAT FORCEPS

12.5cm Straight | SKU 1221



PRODUCT OVERVIEW

Summary	1 x 2 teeth. Straight. Overall length 5…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Kocher Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 12.5 cm Tip Type: 1x2 Teeth & Horizontally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Kocher's artery forceps, kocker's artery forceps, Rochester-Ochsner hemostatic tissue forceps, Ochsner forceps Surgical Instrument: Hospitime kocher's hemostat forceps is a ratcheted forceps utilized to grasp medium to heavy tissue & dense vessels. The instrument is horizontally serrated on the entire length of the jaw . The combination of serrations in addition to teeth guarantee a firm hold on the tissue or veins being held.

PRODUCT TAGS

Kocher's Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips