KOCHER'S HEMOSTAT FORCEPS

15cm Curved | SKU 1224



PRODUCT OVERVIEW

Summary	1 x 2 teeth. Curved. Overall length 6…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Kocher Instrument Type: Hemostatic Forceps Blade Type: Curved Length: 15 cm Tip Type: 1x2 Teeth & Horizontally Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Kocher's artery forceps, kocker's artery forceps, Rochester-Ochsner tissue forceps, Ochsner forceps Surgical Instrument: Hospitime kocher's hemostat forceps is a ratcheted forceps utilized to grasp medium to heavy tissue & dense vessels. The instrument is horizontally serrated on the entire length of the jaw . The combination of serrations in addition to teeth guarantee a firm hold on the tissue or veins being held.

PRODUCT TAGS

Kocher's Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips