LANE'S TWIN INTESTINAL CLAMP STRAIGHT

SKU 1338

[product image not available]

PRODUCT OVERVIEW

Summary	Longitudnal serration. Straight.
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Instrument Type: Hemostatic Clamps Blade Type: Straight Tip Type: Longitudinally Serrated Handle: 2 separate handles with Ratchet & Finger Rings Material: Stainless Steel Alias Names: Lane twin stomach clamp, Lanes twin gastrojejunostomy intestinal occlusion clamp, Tanner twin clamp Surgical Instrument: Hospitime lane's twin intestinal clamp is used for gastrojejunal anastomosis. It has two pairs of blades one each for stomach and small intestine and are non-crushing type.

PRODUCT TAGS

Lane's Twin Intestinal Clamp Straight | General Surgery Instruments | Forceps, Clamps & Clips