LANGENBECK RETRACTOR

20mm Width | SKU 1430

[product image not available]

PRODUCT OVERVIEW

Summary	Right angle blade, blade length 20m…
Category	General Surgery Instruments, Retra…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Langenbeck Instrument Type: Retractors Blade Width: 20X60mm Blade Type: Right Angled Flat Blade Length: 22 cm Tip Type: Curved Ends Handle: Ring Handle Material: Stainless Steel Alias Names: Right angle retractor Surgical Instrument: Hospitime langenbeck is a basic handheld instrument with a flat blade bent down at right angles to the handle. This surgical instrument is helpful for retraction of soft and delicate tissues during the procedures.

PRODUCT TAGS

Langenbeck Retractor | General Surgery Instruments | Retractor