
LIGATURE HEMOSTAT FORCEPS "MIXTER"

23cm | SKU 1248



PRODUCT OVERVIEW

Summary	Fully curved. Overall length 9in (23…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Mixter Instrument Type: Hemostatic Forceps Blade Type: Right Angled (90Â° Degrees) Length: 23 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Ligature artery forceps, Mixter artery forceps, Ligature forceps, right angle artery forceps Surgical Instrument: Hospitime ligature hemostat forceps (mixter) is a ratcheted forceps utilized for working in obscured surgical locales. They are generally utilized for clamping, dissection & handling tissue. Mixter right angle forceps have serrations on the entire length of the jaw.

PRODUCT TAGS

Ligature Hemostat Forceps "Mixter" | General Surgery Instruments | Forceps, Clamps & Clips