MAYO HEGAR NEEDLE HOLDER

12.5cm | SKU 1352

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 5in (12.5cm).
Category	General Surgery Instruments, Needl…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Mayo-Hegar Length: 12.5 cm Instrument Type: Needle Holder Curvature: Straight Blade Type: Cross-serrated with groove Tip Type: Blunt Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Mayo needle holder Surgical Instrument: Hospitime mayo hegar needle holder is a ratcheted instrument, used in multiple types of surgical techniques to hold a suturing needle for closing wounds during suturing and surgical procedures.

PRODUCT TAGS

Mayo Hegar Needle Holder | General Surgery Instruments | Needle Holder