# MAYO HEGAR NEEDLE HOLDER

20cm | SKU 1355

[ product image not available ]

## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Overall length 8in (20cm). |
| **Category** | General Surgery Instruments, Needl… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Mayo-Hegar Length: 20 cm Instrument Type: Needle Holder Curvature: Straight Blade Type: Cross-serrated with groove Tip Type: Blunt Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Mayo needle holder Surgical Instrument: Hospitime mayo hegar needle holder is a ratcheted instrument, used in multiple types of surgical techniques to hold a suturing needle for closing wounds during suturing and surgical procedures.

## PRODUCT TAGS

Mayo Hegar Needle Holder | General Surgery Instruments | Needle Holder