MAYO'S SCISSORS

16.5cm Curved | SKU 1055



PRODUCT OVERVIEW

Summary	Curved. Overall length 6-1/2in (16.5c…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Mayo Instrument Type: Scissors Blade Type: Curved Length: 16.5 cm Tip Type: Semi-Blunt Ends Handle: Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime mayo scissors are designed for cutting body tissues near the surface of a wound, allows deeper penetration into the wound, used to cut thick tissues and also used for cutting sutures.

PRODUCT TAGS

Mayo's Scissors | General Surgery Instruments | Scissors