MAYO'S SCISSORS

21.5cm Straight | SKU 1058



PRODUCT OVERVIEW

Summary	Straight. Overall length 8-1/2in (21c…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Mayo Instrument Type: Scissors Blade Type: Straight Length: 21.5 cm Tip Type: Semi-Blunt Ends Handle: Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime mayo scissors are designed for cutting body tissues near the surface of a wound, allows deeper penetration into the wound, used to cut thick tissues and also used for cutting sutures.

PRODUCT TAGS

Mayo's Scissors | General Surgery Instruments | Scissors