MC-INDOE DISSECTING FORCEPS

15cm 1X2 Teeth | SKU 1151

[product image not available]

PRODUCT OVERVIEW

Summary	Fine, 1 x 2 teeth. Overall length 6in (…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: McIndoe Instrument Type: Tissue Forceps Curvature: Straight Length: 15 cm Tip Type: 1x2 Teeth (Tooth) Handle: Thumb Material: Stainless Steel Alias Names: McIndoe tissue forceps Surgical Instrument: Hospitime Mc-Indoe forceps 15cm (tooth) are utilized to grasp tissues in a surgical wound, to remove debris and to place or remove dressings & packing. They are thumb forceps with a stop peg to guarantee that when utilized for grasping tissue the grip is not crushing.

PRODUCT TAGS

Mc-Indoe Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Plastic / Dermatology Instruments | Forceps & Clamps