METZENBAUM SCISSORS

17.5cm Curved | SKU 1068



PRODUCT OVERVIEW

Summary	Blunt/Blunt. Curved. Overall length 7…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Metzenbaum Instrument Type: Scissors Blade Type: Curved Length: 17.5 cm Tip Type: Semi-Blunt Ends Handle: Finger Rings Material: Stainless Steel Alias Names: Tonsil scissor Surgical Instrument: Hospitime metzenbaum scissors are designed for cutting delicate tissue and blunt dissection and are not recommended for cutting heavy tissues or sutures. The scissors come in variable lengths and have a relatively long shank-to-blade ratio.

PRODUCT TAGS

Metzenbaum Scissors | General Surgery Instruments | Scissors | ENT Instruments | Tonsillectomy