# METZENBAUM SCISSORS

20cm Curved | SKU 1070




## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | Blunt/Blunt. Curved. Overall length 8… |
| **Category** | General Surgery Instruments, Scisso… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Pattern: Metzenbaum Instrument Type: Scissors Blade Type: Curved Length: 20 cm Tip Type: Semi-Blunt Ends Handle: Finger Rings Material: Stainless Steel Alias Names: Tonsil scissor Surgical Instrument: Hospitime metzenbaum scissors are designed for cutting delicate tissue and blunt dissection and are not recommended for cutting heavy tissues or sutures. The scissors come in variable lengths and have a relatively long shank-to-blade ratio.

## PRODUCT TAGS

Metzenbaum Scissors | General Surgery Instruments | Scissors | ENT Instruments | Tonsillectomy