METZENBAUM SCISSORS

22.5cm Straight | SKU 1071



PRODUCT OVERVIEW

Summary	Blunt/Blunt. Straight. Overall length…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Metzenbaum Instrument Type: Scissors Blade Type: Straight Length: 22.5 cm Tip Type: Semi-Blunt Ends Handle: Finger Rings Material: Stainless Steel Alias Names: Tonsil scissor, Nelson scissor Surgical Instrument: Hospitime metzenbaum scissors are designed for cutting delicate tissue and blunt dissection and are not recommended for cutting heavy tissues or sutures. The scissors come in variable lengths and have a relatively long shank-to-blade ratio.

PRODUCT TAGS

Metzenbaum Scissors | General Surgery Instruments | Scissors