
MICRO SPRING NEEDLE HOLDER WITHOUT LOCK

15cm Straight | SKU 1377



PRODUCT OVERVIEW

Summary	Straight. Without lock. Overall length…
Category	Ophthalmic Instruments, Needle Hol…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Castroviejo Length: 15 cm Instrument Type: Needle Holder Blade Type: Straight Tip Type: Blunt Handle: Thumb (double spring) Material: Stainless Steel Alias Names: Castroviejo microsurgical needle holder Surgical Instrument: Hospitime micro spring needle holder is a double spring instrument utilized for holding small needles in microsurgical procedures. It also features a locking mechanism.

PRODUCT TAGS

Micro Spring Needle Holder Without Lock | Ophthalmic Instruments | Needle Holders | Micro Surgery Instruments | Needle Holder