
MOSQUITO HEMOSTAT FORCEPS

7.5cm Straight | SKU 1170



PRODUCT OVERVIEW

Summary	Straight. Overall length 3in (7.5cm).
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Mosquito Instrument Type: Hemostatic Forceps Blade Type: Straight Length: 7.5 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Mosquito artery forceps, Halsted mosquito forceps, Hartmann mosquito forceps Surgical Instrument: Hospitime mosquito hemostat forceps 7.5cm (straight) are ratcheted forceps which are utilized in different surgical procedures as a hemostat for small blood veins. They are not appropriate for use with large vessels or arteries. They are utilized for holding delicate tissues in skin grafting, biopsies or ophthalmologic methodology. Mosquito forceps are known for their fine tips, and short, completely serrated jaws.

PRODUCT TAGS

Mosquito Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips | Ophthalmic Instruments | Forceps & Clamps