MOYNIHAN HEMOSTAT FORCEPS

SKU 1238



PRODUCT OVERVIEW

Summary	Angled.
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Moynihan Instrument Type: Hemostatic Forceps Blade Type: Curved Tip Type: Blunt Handle: Ratchet & Finger Rings Material: Stainless Steel Surgical Instrument: Hospitime moynihan hemostat forceps are ratcheted utilized for handling tissue or passing sutures in abdominal procedures. They have long shanks and short jaws with horizontal serrations.

PRODUCT TAGS

Moynihan Hemostat Forceps | General Surgery Instruments | Forceps, Clamps & Clips