NEGUS HEMOSTAT FORCEPS

19cm | SKU 1257

[product image not available]

PRODUCT OVERVIEW

Summary	Overall length 7-1/2in (19cm).
Category	General Surgery Instruments
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Negus Instrument Type: Hemostatic Forceps Blade Type: Small Hook Shaped Length: 19 cm Tip Type: Serrated Handle: Ratchet & Finger Rings Material: Stainless Steel Alias Names: Negus artery forceps, Clamp forceps, Tonsil forceps, Negus tonsil artery forceps Surgical Instrument: Hospitime Negus Hemostatic Forceps are ratcheted and finger ring forceps. This instrument replaces tonsil artery forceps; utilized to tie ligature at a depth and ligature would not slip because of its small hook shaped tip . It is curved with an overall length of 7Â½ inches.

PRODUCT TAGS

Negus Hemostat Forceps | General Surgery Instruments