
# PAYR'S CRUSHING INTESTINAL CLAMP BLADE

15 cm | SKU 1333




## PRODUCT OVERVIEW

| | |
|---|---|
| **Summary** | With pin. Blade length 15cm. |
| **Category** | General Surgery Instruments, Force… |
| **Availability** | In stock |

## PRODUCT DESCRIPTION

Instrument Type: Hemostatic Clamps Curvature: Angled Blade Length: 15 cm Tip Type: Longitudinally Serrated Handle: Finger Rings Material: Stainless Steel Alias Names: Payr pylorus clamp, Crushing clamp Surgical Instrument: Hospitime payr's crushing intestinal clamp is utilized to occlude the lumen of the intestine to avoid spillage of gut content and to occlude the blood supply to the diseased part of the gut. The forceps are available with and without a pin depending on doctor's demand. The clamp has four joints, the three joints in the vicinity are for the handles and they operate the distal joint. The instrument is double action with an assortment of jaw lengths.

## PRODUCT TAGS

Payr's Crushing Intestinal Clamp Blade | General Surgery Instruments | Forceps, Clamps & Clips