
POTT'S SCISSORS 60° ANGLED

SKU 1093



PRODUCT OVERVIEW

Summary	Angled 60°. Overall length 7in (17.5…
Category	General Surgery Instruments, Scisso…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Pott's Instrument Type: Scissors Blade Type: 60° Angled Length: 17.8 cm Tip Type: Sharp Handle: Finger Rings Material: Stainless Steel Alias Names: Angled scissors, Potts-Smith scissors Surgical Instrument: Hospitime pott's scissors is an angled scissor which is significantly utilized for trimming and opening vessels in cardiovascular and thoracic methodology.

PRODUCT TAGS

Pott's Scissors 60° Angled | General Surgery Instruments | Scissors | Cardio Thoracic & Vascular