RUSSIAN DISSECTING FORCEPS

15cm | SKU 1124



PRODUCT OVERVIEW

Summary	Russian style tip. Overall length 6in (…
Category	General Surgery Instruments, Force…
Availability	In stock

PRODUCT DESCRIPTION

Pattern: Russian Instrument Type: Tissue Forceps Curvature: Straight Length: 15 cm Tip Type: Cup-shaped and radially serrated Handle: Thumb Material: Stainless Steel Alias Names: Russ Modell forceps, Rochester tissue forceps, Russian tissue forceps Surgical Instrument: Hospitime russian Forceps are utilized for handling heavy or thick tissue. They are likewise utilized during wound closure. This item has atraumatic, radially serrated and cupped tips.

PRODUCT TAGS

Russian Dissecting Forceps | General Surgery Instruments | Forceps, Clamps & Clips