SCALPEL HANDLE B.P NO. 3L

SKU 1003



PRODUCT OVERVIEW

Summary	For blades nos. 10, 11, 12 and 15.
Category	General Surgery Instruments, Scal…
Availability	In stock

PRODUCT DESCRIPTION

Uses size 10 - 15 blades Material: Stainless Steel Alias Names: Surgical knife, BardParker handle, BP handle, Dissection knife Surgical Instrument: Knife handles are used to hold various blades to create scalpel. Scalpels are used to make skin incisions or whenever a fine precision cut is necessary. Hospitime Scalpel Handle #3L is used for surgery & anatomical dissection. It is longer then Scalpel handle #3.

PRODUCT TAGS

Scalpel Handle B.P No. 3L | General Surgery Instruments | Scalpel, Blade breaker & Caliper