SCALPEL HANDLE B.P NO. 4L

SKU 1004-1



PRODUCT OVERVIEW

Summary	For blades nos. 20, 21, 22, 23, 24 a…
Category	General Surgery Instruments, Scal…
Availability	In stock

PRODUCT DESCRIPTION

Uses size 20 - 25 blades. Material: Stainless Steel Alias Names: Surgical knife, BardParker handle, BP handle, Dissection knife Surgical Instrument: Knife handles are used to hold various blades to create scalpel. Scalpels are used to make skin incisions or whenever a fine precision cut is necessary. Hospitime Scalpel Handle #4L is used for surgery & anatomical dissection. It is longer then Scalpel handle #4.

PRODUCT TAGS

Scalpel Handle B.P No. 4L | General Surgery Instruments | Scalpel, Blade breaker & Caliper